

# Sports Medicine product catalog



[actimovesport.com](http://actimovesport.com)

**Actimove**

**Leukoplast**

**Delta-Cast**

**JOBST**

 **essity**

# Forward

It's just not a direction.  
**It's an attitude.**

It's about not resting on our laurels,  
**but leading the way.**

**Forward-thinking is what makes us a leader in health, hygiene, and medical solutions.**

Our extensive portfolio of leading therapeutic brands reflects the spirit of innovation we bring to everything we do.

We have the courage to lead, to break barriers, and to always push for better.

**Better Products.**

**Better Service.**

**Better Value.**

**Better Outcomes.**

Our wide range of therapeutic solutions, clinical expertise and world-class brands **mean one partner to support your needs.**

Orthopedics: **Actimove** **Delta-Cast**

Wound Care: **Leukoplast**

Compression: **JOBST**

## Your advantage

### Convenience

- One partner for a range of therapeutic solutions
- The right assortment tailored to your needs

### Medical expertise

- On-going educational tools and support
- Quality products based on validated evidence for better outcomes
- Confidently recommend the right product

## Our goal

**Being your ally in better outcomes for you and those in your care.  
Now and for the future.**



# Contents

## Sport Taping & Bandaging

Introduction	6
Portfolio overview	7
<b>1.1 ELASTIC ADHESIVE TAPES</b>	
ELASTIKON® by Actimove®	8
TensoSport®	8
Tensoplast®	9
Tensoplast® LF	9
<b>1.2 LIGHT ELASTIC ADHESIVE BANDAGES</b>	
Lightplast® Pro	10
<b>1.3 NON-ELASTIC ADHESIVE STRAPPING TAPES</b>	
Leukotape® P	10
COACH® by Actimove®	11
ZONAS® by Actimove®	11
<b>1.4 ADHESIVE PROTECTION SPRAY</b>	
Tensospray®	11
<b>1.5 KINESIOLOGY TAPES</b>	
Leukotape® K	12
<b>1.6 NON-ADHESIVE SUPPORT BANDAGES</b>	
BSN Elastic Bandage	12
<b>1.7 COHESIVE LIGHT SUPPORT BANDAGES</b>	
Co-Plus® LF	13
<b>1.8 TUBULAR SUPPORT BANDAGES</b>	
Tensogrip®	13

## Fracture Management

Introduction	14
Portfolio overview	15
<b>2.1 CAST TAPE</b>	
Delta-Lite® Plus	16
Delta-Lite® Conformable	16
<b>2.2 PADDINGS</b>	
Delta-Rol® Synthetic Cast	17
Specialist® Micro Pleat Cotton Blend Cast	17
Sof-Rol® Absorbent Rayon Cast	18
Delta-Dry® Water Resistant Padding	18
<b>2.3 STOCKINETTE</b>	
AquaCast® Liner	19
Delta-Net® Orthopedic Synthetic	20
Specialist® Orthopedic Cotton	20
Terry-Net™ S Synthetic Terry Cloth	20
Terry-Net™ Thumb Spica	21
<b>2.4 CAST ACCESSORIES</b>	
Delta-Net® Adhesive Felt Padding	21
Sifoam™ Interface Padding	21

## Fracture Management (cont.)

Delta Terry-Net™ Adhesive Fleece Edger	22
Delta Terry-Net™ Adhesive Hook & Stretch Loop Closures	22
Zip Stick	22
American Orthopaedic™	23
Essity Buckets	23
Saw Stop® DE-FLEX™	23
Stainless Steel Cast Spreader	24
FCT Scissors	24
Esmarch Scissors	24
ORTHO-GLASS® Scissors	24
<b>2.5 SPLINTING</b>	
ORTHO-GLASS®	25

## Orthopedic Soft Goods

Introduction	26
Portfolio overview	27
<b>3.1 Actimove® UPPER EXTREMITY SUPPORTS</b>	
Actimove® Shoulder Support	28
Actimove® Mitella Comfort Arm Sling	28
Actimove® Elbow Strap	29
Actimove® Elbow Support	29
Actimove® Episport® Elbow Clasp	29
Adjustable	29
Actimove® Wrist Stabilizer	30
Actimove® Thumb Stabilizer	30
Actimove® Alufoam	31
Actimove® Rhizo Forte Thumb	31
<b>3.2 Actimove® LOWER EXTREMITY SUPPORTS</b>	
Actimove® Knee Brace	32
Actimove® Knee Brace	32
Actimove® Knee Brace	33
Actimove® Knee Stabilizer	33
Actimove® Knee Support	34
Actimove® Post Op Rom	34
Actimove® Actimove® Thigh Support	35
Actimove® Thigh Support	36
Actimove® PowerMotion Thigh Support	36
Actimove® PowerMotion Calf Support	37
Actimove® Standard Walker Low / High	37
Actimove® Standard Air Walker Low / High	38
Actimove® Closed Shell Air Walker Low / High	38
Actimove® Talocast Ankle Brace	39
Actimove® Ankle Stabilizer	39
Actimove® Patella Strap	40

## Wound Care

Introduction	42
Portfolio overview	43
<b>4.1 STRETCHABLE LIGHTWEIGHT FABRIC DRESSING</b>	
Leukoplast® elastic	44
<b>4.2 TRANSPARENT DRESSING WITH ABSORBENT PAD</b>	
Leukomed® T Plus	45
<b>4.3 WIDE AREA FIXATION</b>	
Leukoplast® Hypafix®	45
<b>4.4 ELASTIC GAUZE BANDAGE</b>	
Elastomull®	46
<b>4.5 HYDROPHOBIC MICROBE BINDING FILM DRESSING</b>	
Leukomed® Sorbact®	46
<b>4.6 CONFORMING, ADHESIVE FIXATION DRESSING</b>	
Cover-Roll® Stretch	47
<b>4.7 WOUND &amp; EYE IRRIGATION</b>	
SplashCap®	47
EyeCap®	47





# Sport Taping & Bandaging

Essity has a full line of quality sports tapes to help support your athletes.

## Portfolio overview

ELASTIC ADHESIVE TAPES	<p>ELASTIKON® by Actimove®</p>	<p>TensoSport® Tensoplast®</p>	<p>Tensoplast®</p>	<p>Tensoplast® LF</p>	
	LIGHT ELASTIC ADHESIVE BANDAGES	<p>Lightplast® Pro</p>	NON-ELASTIC ADHESIVE STRAPPING TAPES		
		<p>Leukotape® P</p>			<p>COACH® by Actimove®</p>
	ADHESIVE PROTECTION SPRAY	<p>Tensospray®</p>	KINESIOLOGY TAPES	NON-ADHESIVE SUPPORT BANDAGES	
<p>Leukotape® K</p>		<p>BSN Elastic Bandage</p>			
COHESIVE LIGHT SUPPORT BANDAGES	<p>Co-Plus® LF</p>			TUBULAR SUPPORT BANDAGES	
				<p>Tensogrip®</p>	

## ELASTIKON® by Actimove®

Soft cotton elastic tape, centerline for ease of application



- Reliable compression for support of sprains, strains and muscle injuries
- Soft cotton elastic tape conforms to joints such as elbows, knees and shoulders
- Allows skin to breathe and moisture to pass through
- Treatment of conditions which require limitation of movement, e.g. sprains, sports injuries
- Prophylactic use in active sport
- Retention of dressings and devices, e.g. post-orthopedic and general surgery or over a muscular wound

Order Number	Size	Color	UOM
5170	2 in. x 2.5 yds	Cream	24 rolls/case (Speed Pack)
5171	3 in. x 2.5 yds	Cream	16 rolls/case (Speed Pack)
5172	1 in. x 2.5 yds	Cream	144 rolls/case (12 rolls/box, 12 box/case)
5174	2 in. x 2.5 yds	Cream	72 rolls/case (6 rolls/box, 12 box/case)
5175	3 in. x 2.5 yds	Cream	48 rolls/case (4 rolls/box, 12 box/case)
5177	4 in. x 2.5 yds	Cream	36 rolls/case (6 rolls/box, 6 box/case)

This Product Contains Dry Natural Rubber.

## TensoSport®

Soft cotton elastic tape, centerline for ease of application

TensoSport® is an elastic adhesive bandage suitable for medium support or compression. Excellent stretch and regain provides firm, controlled compression and support. Fluffy non-fray edges prevent the bandage from rucking or marking.

- High adhesive properties
- Provides strong compression and long-lasting support
- Water repellent

### TensoSport Elastic Adhesive Bandage – Boxed Packs

Order Number	Size	Color	UOM
7205015	2 in. x 5 yds	White	32 rolls/case
7205016	3 in. x 5 yds	White	24 rolls/case

Caution: This product contains natural rubber latex which may cause allergic reactions



## Tensoplast®

Soft cotton elastic tape

Tensoplast® is an elastic adhesive bandage suitable for medium support or compression. Excellent stretch and regain provides firm, controlled compression and support. Fluffy non-fray edges prevent the bandage from rucking or marking.

- High adhesive properties
- Provides strong compression and long-lasting support
- Water repellent



### Tensoplast® Elastic Adhesive Bandage – Individually Boxed

Order Number	Size	Color	UOM
02598002	1 in. x 5 yds	Beige	36 rolls/case
02599002	2 in. x 5 yds	Beige	36 rolls/case
02600002	3 in. x 5 yds	Beige	36 rolls/case
02601002	4 in. x 5 yds	Beige	36 rolls/case
02602002	6 in. x 5 yds	Beige	12 rolls/case
02593002	1 in. x 5 yds	White	36 rolls/case
02594002	2 in. x 5 yds	White	36 rolls/case
02595002	3 in. x 5 yds	White	36 rolls/case
02596002	4 in. x 5 yds	White	36 rolls/case
02597002	6 in. x 5 yds	White	12 rolls/case

### Tensoplast® Elastic Adhesive Bandage – Bulk Packs

Order Number	Size	Color	UOM
04410001	1 in. x 5 yds	Beige	32 rolls/case
04412001	2 in. x 5 yds	Beige	24 rolls/case
04413001	3 in. x 5 yds	Beige	16 rolls/case
04414001	4 in. x 5 yds	Beige	16 rolls/case

Caution: This product contains natural rubber latex which may cause allergic reactions

## Tensoplast® LF

Tensoplast® LF is not made with natural rubber latex and assures the quality you trust. It provides flexible, strong support and compression thanks to its porously high-quality substrate. Its excellent stretch and adhesive molds and holds to any part of the body, keeping the bandages firmly in place even throughout lengthy treatment. Trusted patient comfort and compliance for those sensitive to latex make Tensoplast® LF the optimal choice for a wide range of indications.

- Permeable material
- High-quality elastic substrate
- Adhesive coating
- Yellow center line
- Not made with natural rubber latex

Order Number	Size	Color	UOM
7623200	2 in. x 4.9 yds	Cream	12 rolls/case
7623201	3 in. x 4.9 yds	Cream	12 rolls/case
7623202	4 in. x 4.9 yds	Cream	12 rolls/case



1.2 LIGHT ELASTIC ADHESIVE BANADAGES

## Lightplast® Pro

### Tearable Adhesive Stretch Bandage

Lightplast® Pro is a tearable, elastic adhesive bandage, suitable for light support to compression. The unique fabric backing helps to conform to body movement. Ideal for all-purpose taping.

- Hand-tearable
- Excellent adhesion
- Light and close fitting

Order Number	Size	Color	UOM
76952	1 in. x 5 yds	White	48 rolls/case
76953	1.5 in. x 5 yds	White	32 rolls/case
76954	2 in. x 5 yds	White	24 rolls/case
76955	3 in. x 5 yds	White	16 rolls/case
76956	2 in. x 7.5 yds	White	24 rolls/case
76957	3 in. x 7.5 yds	White	16 rolls/case
72880	2 in. x 5 yds	Black	24 rolls/case
72881	3 in. x 5 yds	Black	16 rolls/case

This Product Contains Dry Natural Rubber.



1.3 NON-ELASTIC ADHESIVE STRAPPING TAPES

## Leukotape® P

### Rigid Strapping Tape

Leukotape® P is a high strength rigid strapping tape with a very strong zinc oxide adhesive. It is porous yet offers strong support for sprains and strains.

- Porous
- Easy unwinding tension
- High strength adhesive
- Ideal for patellofemoral taping

Order Number	Measurements	Color	UOM
76168	1.5 in. x 15 yds	Beige	30 rolls/case

This Product Contains Dry Natural Rubber.



## COACH® by Actimove®

### Rigid Cloth Porous Adhesive Strapping Tape



- Support and conformity helps protect against sprains and injury
- High tensile strength
- Adheres securely yet tears easily by hand
- 100% cotton back cloth allows perspiration to escape
- Smooth and consistent unwind

Order Number	Measurements	Color	UOM
5188	1.5 in. x 15 yds	White	32 rolls/case (Speed Pack)

This Product Contains Dry Natural Rubber.

## ZONAS® by Actimove®

### Rigid Cloth Porous Adhesive Strapping Tape



- Prophylactic use in active sport
- Ideal for conditions requiring immobilization and support
- Highest tensile strength
- Adheres securely yet tears by hand
- 100% cotton back cloth allows perspiration to escape
- Smooth and consistent unwind

Order Number	Measurements	Color	UOM
5104	1 in. x 10 yds	White	144 rolls/case (12 rolls/box, 12 box/case)
5190	1.5 in. x 15 yds	White	32 rolls/case (Speed Pack)

This Product Contains Dry Natural Rubber.

1.4 ADHESIVE PROTECTION SPRAY

## Tensospray®

### Adhesive Protection Spray

Tensospray® is an adhesive spray which leaves a thin protective film over the surface of the skin. This provides an ideal layer for adhesive bandages and tapes, promoting excellent traction and reducing skin trauma on removal. The product is not affected by disinfectant.

- Aerosol presentation makes for quick and easy application
- Provides a copolymer layer which improves the effectiveness of the strapping
- Protects skin and encloses hairs – aids comfort and reduces trauma of bandage removal

Order Number	Measurements	UOM
7160200	300ml	12 box





1.5 KINESIOLOGY TAPES

# Leukotape® K

## Kinesiology Therapeutic Adhesive Tape

The wave pattern adhesive of **Leukotape® K** allows air to circulate. It is stretchable and works ideally with the skin's natural elasticity.

- High wearing comfort and skin friendly, with a formulation not made with natural rubber latex
- Stays in place up to 9 days, even after showering or sweating
- Residue free removal
- Dermatologically tested



Order Number	Size	Color	UOM / Pack
7297810	1 in. x 5.5 yds		
7297811	2 in. x 5.5 yds	● Tan	
7297812	3 in. x 5.5 yds		
7297815	1 in. x 5.5 yds		
7297816	2 in. x 5.5 yds	● Red	
7297817	3 in. x 5.5 yds		5 rolls/pack
7297820	1 in. x 5.5 yds		
7297821	2 in. x 5.5 yds	● Blue	
7297822	3 in. x 5.5 yds		
7297823	2 in. x 5.5 yds	● Black	
7297824	2 in. x 5.5 yds	● Lt Blue	
7297825	2 in. x 5.5 yds	● Pink	

1.6 NON-ADHESIVE SUPPORT BANDAGES

# BSN Elastic Bandage

## Non-Adhesive Bandage

- Adjustable and constructed of breathable elastic
- Hook and loop closure
- No metal clips
- Not made with natural rubber latex
- Individually wrapped



Order Number	Size	Color	UOM / Pack
1037023	2 in. x 5 yds		
1037033	3 in. x 5 yds		
1037043	4 in. x 5 yds	● Beige	10 rolls/pack, 5 packs/case
1037053	6 in. x 5 yds		

1.7 COHESIVE LIGHT SUPPORT BANDAGES

# Co-Plus® LF

## Cohesive Bandage

- Self-adhesive
- Non-slip support
- Hand-tearable in both directions
- End-of-bandage indicator for ease of use
- Not made with natural rubber latex
- Available in beige and color pack: neon pink, light blue, purple, neon green and red



Order Number	Size	Color	UOM
7210014	1 in. x 5 yds		60 rolls/box
7210016	2 in. x 5 yds		36 rolls/pack
7210017	3 in. x 5 yds	● Beige	24 rolls/pack
7210018	4 in. x 5 yds		18 rolls/pack
7210019	6 in. x 5 yds		12 rolls/pack
7210020	3 in. x 5 yds		24 rolls/pack
7210021	4 in. x 5 yds	Assorted	18 rolls/pack
7210023	2 in. x 5 yds		36 rolls/pack

1.8 TUBULAR SUPPORT BANDAGES

# Tensogrip®

## Easy to apply and help control edema

- Ideal for the light support of strains, sprains and soft tissue injuries
- Provides light support, light compression and dressing retention
- Constructed from a cotton/polyamide blend with natural rubber latex yarn



Order Number	Size	Color	UOM
7579	A - 1.75 in x 11 yds	White	1/roll
7580	B - 2.5 in x 11 yds	White	1/roll
7581	C - 3 in x 11 yds	White	1/roll
7582	D - 3.5 in x 11 yds	White	1/roll
7583	E - 3.5 in x 11 yds	White	1/roll
7584	F - 4 in x 11 yds	White	1/roll
7585	G - 4.5 in x 11 yds	White	1/roll
7586	J - 7 in x 11 yds	White	1/roll
7587	K - 8.5 in x 11 yds	White	1/roll

Caution: This product contains natural rubber latex which may cause allergic reactions



# Portfolio overview

**CSAT TAPE**

Delta-Lite® Plus    Delta-Lite® Conformable

**PADDINGS**

Delta-Rol®    Specialist®    Sof-Rol®    Delta-Dry®    AquaCast® Liner

**STOCKINETTE**

Delta-Net®    Specialist®    Terry-Net™    Terry-Net™ Thumb Spica

**CAST ACCESSORIES**

Delta-Net®    Sifoam™    Delta Terry-Net™    Delta Terry-Net™    Zip Stick

American Orthopaedic™    Essity® Buckets    Saw Stop® DE-FLEX™    Stainless Steel Cast Spreader    FCT Scissors

Esmarch Scissors    ORTHO-GLASS® Scissors

**SPLINTING**

ORTHO-GLASS®

## Fracture Management

Essity has a full line of casting and splinting material to support your orthopedic injuries.



## Delta-Lite® Plus

### The Trusted Rigid Solution Fiberglass Cast Tape

**Delta-Lite® Plus** is designed to ensure excellent rigid support. Thanks to its light but strong fiberglass substrate and resin formula it has become a trusted solution for fractures requiring strong and lasting support.

- Low tack resin
- Primary and secondary casting
- Rigid support
- Good conformability
- High lamination
- Excellent end lay down
- Smooth surface
- 3 – 5 minute initial set time
- Available in 11 colors and mixed packs
- Weight bearing in 20 minutes

Order Number	Size	Color	UOM
7345800	1 in x 4 yds	White	10 rolls/box
7345801	2 in x 4 yds	White	10 rolls/box
7345802	3 in x 4 yds	White	10 rolls/box
7345803	4 in x 4 yds	White	10 rolls/box
7345804	5 in x 4 yds	White	10 rolls/box



## Delta-Cast® Conformable

### The Non-Fiberglass Solution in Cast Tape

Sets a benchmark in non-fiberglass cast tapes. **Delta-Cast® Conformable** responds to all challenges of daily fracture management. The secret behind Delta-Cast® Conformable and its outstanding benefits lies in its unique components, making it the premium choice in cast rooms around the world. Provides rigid and semi-rigid cast-ing options for primary and secondary casting applications where patient fit, comfort and compliance are critical.

- Primary and secondary casting
- Rigid and semi-rigid support
- Excellent conformability
- Weight bearing in 20 minutes
- 3 – 5 minute initial set time
- 9 colors and mixed packs
- Excellent for Functional Cast Therapy
- Excellent radiolucency
- Fiberglass-free
- Use the gloves of your choice
- Strong and lightweight
- Soft edges
- Excellent lamination
- Smooth finish

Order Number	Size	Color	UOM
6601	1 in x 2 yds	White	10 rolls/box
6602	2 in x 4 yds	White	10 rolls/box
6603	3 in x 4 yds	White	10 rolls/box
6604	4 in x 4 yds	White	10 rolls/box
6605	5 in x 4 yds	White	10 rolls/box



## Delta-Rol®

### Synthetic Cast Padding

**Delta-Rol® Synthetic Cast Padding** is a non-absorbent synthetic material that does not hold moisture against the skin. Its conformable stretch allows narrow widths around small anatomies without cutting and tearing.

- Conformable stretch eliminates need for scissors
- Comfortable
- Non-absorbent synthetic fiber

Order Number	Size	Color	UOM
6882	2 in x 4 yds	White	72 rolls/box
6883	3 in x 4 yds	White	36 rolls/box
6884	4 in x 4 yds	White	36 rolls/box
6886	6 in x 4 yds	White	18 rolls/box



## Specialist®

### Micro Pleat Cotton Blend Cast Padding

**Specialist® Micro Pleat Cotton Blend Cast Padding** has a unique micropleated design for excellent conformability, protection and ease of application without tucking or tearing.

- Micropleated construction provides excellent conformability
- Cotton rayon blend
- Easy application

Order Number	Size	Color	UOM
9062	2 in x 4 yds	White	144 rolls/box
9043	3 in x 4 yds	White	72 rolls/box
9044	4 in x 4 yds	White	72 rolls/box
9046	6 in x 4 yds	White	36 rolls/box



## Sof-Rol®

### Absorbent Rayon Cast Padding

Sof-Rol® is a highly absorbent Rayon material suitable for use around wounds. This soft natural material creates a comfortable undercast surface for either synthetic or natural casting materials.

- Excellent conformability
- Strong lamination
- Excellent molding “touch”

Order Number	Size	Color	UOM
9052	2 in x 4 yds	White	144 rolls/box
9033	3 in x 4 yds	White	72 rolls/box
9034	4 in x 4 yds	White	72 rolls/box
9036	6 in x 4 yds	White	36 rolls/box



## Delta-Dry®

### Water Resistant Padding

The patented Delta-Dry® material offers water resistance. The open knit structure of the fibers accelerate drainage after submerging in water. The patient’s body temperature will then evaporate residual moisture quickly through the airy mesh (high MVTR). In this way, cast and skin dry quickly and skin feels dry in about 90 minutes after daily hygiene.

- Allows patients to shower, swim, bathe and wash hands
- Helps reduce odor under the cast
- Excellent conformability to body contours
- Provides comfortable, fast-drying casts
- Not made with natural rubber latex

Order Number	Size	Color	UOM
7344300	2 in x 2.6 yds	White	12 rolls/box
7344301	3 in x 2.6 yds	White	12 rolls/box
7344302	4 in x 2.6 yds	White	12 rolls/box



## AquaCast® Liner

### Waterproof Cast Padding

AquaCast® Liner is a comfortable waterproof, breathable cast padding for use under waterproof casting materials. It is moisture vapor permeable and made with pillows and channels to enhance comfort and facilitate drainage.

- Quick drying of cast allows patients to bathe, swim and exercise
- Decrease odor and enhance hygiene with the ability to thoroughly wash and rinse skin
- Supports favorable skin condition following cast removal
- Less itching due to smooth material surfaces
- Less recasting due to wet or damaged casts
- Not made with natural rubber latex



Order Number	Size	Color	UOM
<b>Aquacast® Liner Rolls</b>			
ACL-2-S	2 in	White	12 roll/case
ACL-3-S	3 in	White	12 roll/case
ACL-4-S	4 in	White	12 roll/case
ACL-6-S	6 in	White	6 roll/case
<b>Aquacast® Form Fit Liners</b>			
FFS	Form Fit Liner - Small	White	4 units / bag
FFM	Form Fit Liner - Medium	White	4 units / bag
FFL	Form Fit Liner - Large	White	4 units / bag
<b>Aquacast® Liner Starter Pack and Cast Kits</b>			
ACLSP	Starter Kit	White	1 kit
2ACL-2-S	Cast Kit 2": 2 rolls 2", 1 roll SawStop	White	1 kit
3ACL-2-S	Cast Kit 2": 3 rolls 2", 1 roll SawStop	White	1 kit
2ACL-3-S	Cast Kit 3": 2 rolls 3", 1 roll SawStop	White	1 kit
3ACL-3-S	Cast Kit 3": 3 rolls 3", 1 roll SawStop	White	1 kit
2ACL-4-S	Cast Kit 4": 2 rolls 4", 1 roll SawStop	White	1 kit
3ACL-4-S	Cast Kit 4": 3 rolls 4", 1 roll SawStop	White	1 kit

The AquaCast® Liner Starter Pack includes 2 rolls of 2", 3", 4" AquaCast® Liner plus a full-length roll of SawStop® DE-FLEX® Protective Strip.

## Delta-Net®

### Orthopedic Synthetic Stockinette

- Soft and conformable
- Used for Casting and under bracing

Order Number	Size	Color	UOM
6861	1 in x 25yds	White	2 rolls/case
6862	2 in x 25yds	White	2 rolls/case
6863	3 in x 25yds	White	2 rolls/case
6864	4 in x 25yds	White	2 rolls/case
6866	6 in x 25yds	White	1 rolls/case
6868	8 in x 25yds	White	1 rolls/case
6810	10 in x 25yds	White	1 rolls/case



## Specialist®

### Orthopedic Cotton Stockinette

- Soft and conformable
- Used for Casting and under bracing

Order Number	Size	Color	UOM
9071	1 in x 25yds	White	1 rolls/case
9072	2 in x 25yds	White	9 rolls/case
9073	3 in x 25yds	White	8 rolls/case
9074	4 in x 25yds	White	8 rolls/case
9076	6 in x 25yds	White	6 rolls/case



## Terry-Net™ S

### Synthetic Terry Cloth Stockinette

- Soft and conformable
- Terry-net cushion
- High stretch
- Used for Casting and under bracing

Order Number	Size	Color	UOM
53012	2 in x 15yds	White	2 rolls/case
53013	3 in x 15yds	White	2 rolls/case
53014	4 in x 15yds	White	1 rolls/case
53015	5 in x 15yds	White	1 rolls/case
53016	6 in x 15yds	White	1 rolls/case



## Terry-Net™ Thumb Spica

### Precut Synthetic Terry Cloth Stockinette

- Soft and conformable
- Terry-net cushion
- High stretch
- Used for Casting and under bracing

Order Number	Size	Color	UOM
53120	2 in x 16yds (short arm)	White	10 box/case
53121	2 in x 19yds (long arm)	White	10 box/case
53101	3 in x 16yds (short arm)	White	10 box/case
53102	3 in x 19yds (long arm)	White	10 box/case



## 2.4 CAST ACCESSORIES

### Delta Terry-Net™

#### Adhesive Felt Padding

Delta Terry-Net™ Adhesive Felt Padding is an adhesive orthopedic felt, used for added patient comfort. Its adhesive backing makes attaching the padding simple and easy. It may also be used as a cast edge padding or to provide protection for bony prominences

- Soft orthopedic felt
- Adhesive back

Order Number	Size	Color	UOM
54012	23 in x 39 in	White	1 roll/case



### Sifoam™

#### Interface Padding

- Offloads pressure around bony prominences
- Helps prevent shearing
- Is compressable
- Unloads weight

Order Number	Size	Color	UOM
7279900	8 in x 9 in	Black	10 units/box





## Delta Terry-Net™

### Adhesive Fleece Edger

Delta Terry-Net™ Adhesive Fleece Edger features a high quality adhesive for secure placement. Its unique stretch allows the user to pad easily around corners and contours, ensuring smooth cast edges.

- Highly conformable
- Unique stretch
- Secure adhesive
- Works great with thermoplastics

Order Number	Size	Color	UOM
56012	1.25 in x 15 yds	White	1 roll/case



## Delta Terry-Net™

### Adhesive Hook & Stretch Loop Closures

This easy-to-use product simplifies the FCT technique. It provides closure without creating the discomfort caused by traditional loops' rough edges. Its high quality adhesive supports secure placement and the simple dispensing box provides convenient storage.

- Soft stretch loop
- High quality adhesive
- Dispensing box

Order Number	Size	Color	UOM
<b>Adhesive Hook</b>			
52100	1 in x 10 yds	White	1 roll/case
52200	2 in x 10 yds	White	1 roll/case
<b>Stretch Loop</b>			
52101	1 in x 15 yds	White	1 roll/case
52201	2 in x 15 yds	White	1 roll/case




## Zip Stick

### Cast Removal Aid

This convenient cast aid provides additional patient protection during cast removal. The stick slides under the cast and provides a smooth path for scissors or cast saw. Standard and metric rule adds functionality.

Order Number	Color	UOM
7204629	White	1 ea



## American Orthopaedic™

### Cast Cutter Cast Saws

The American Orthopaedic™ Cast Cutter provides precisely what is needed when needed: power, speed, quiet operation, durability, ease of operation and a balance of weight and power. Includes one 2.5" diameter multi-purpose blade and wrench.

- Balance of weight and power helps reduce user fatigue
- Internally lubricated for longer service
- Designed for rugged use
- Switch is designed and placed so that it will not break if dropped
- Hex drive provides six (6) blade rotations extending blade life by 50%

Order Number	Color	UOM
0295-200	Blue	1/box



## Essity® Buckets

The wax coated bucket is used to hold water for cast activation.

Order Number	Size	Color	UOM
4183-155	5 quart	White	12/box



## Saw Stop® DE-FLEX™

These protective strips are tough, waterproof, flexible ribbons that add an added layer of protection from a cast saw.

- Adhesive Backing for easy application

Order Number	Size	Color	UOM
SS10	1.25 in x 10 ft	Blue	1/box
SS30	1.25 in x 10 ft	Blue	3/box
SS100	1.25 in x 10 ft	Blue	10/box
SSW10	1.25 in x 10 ft	White	1/box
SSW30	1.25 in x 10 ft	White	3/box
SSW100	1.25 in x 10 ft	White	10/box



## Casting Instruments

For the Orthopedic Specialist

Casting instruments from BSN medical are effective for cutting casting material. They are handcrafted with high quality stainless steel to provide precision performance.

### Stainless Steel Cast Spreader

- Spring loaded
- Excellent to open cut casts in cast removal

Order Number	Size	UOM
21-0136	9 in	1 unit



### FCT Scissors

- Sharp, blunt tip scissors
- Clean cut
- Excellent for cast trimming

Order Number	Size	UOM
28230	8 in	1 unit



### Esmarch Scissors

- Sharp, blunt tip scissors
- Excellent to cut bandages

Order Number	Size	UOM
7204407	8 in	1 unit



### ORTHO-GLASS® Scissors

- Sharp, no blunt tip scissors
- Excellent for cutting Splints

Order Number	Size	UOM
7204657	9 in	1 unit



## 2.5 SPLINTING

### ORTHO-GLASS®

Fiberglass Padded Splinting System

ORTHO-GLASS® is a cost-effective splinting system of choice for every phase of immobilization. Features Interlocking Performance Technology (ILP). ILP eliminates the risk of delamination, ensuring a strong, durable splint.

- Eliminates risk of delamination
- No wrong side application
- Easy stretch padding covers all edges
- Roll form eliminates waste
- Full educational support system - workshops, in-services, manuals, wall-charts, videos



#### ORTHO-GLASS® Splinting System

Order Number	Size	UOM
OG-1L2	1 in x 15 ft	2 rolls / case
OG-2L2	2 in x 15 ft	2 rolls / case
OG-3L2	3 in x 15 ft	2 rolls / case
OG-4L2	4 in x 15 ft	2 rolls / case
OG-5L2	5 in x 15 ft	2 rolls / case
OG-6L2	6 in x 15 ft	2 rolls / case
66000687	8 in x 7.5 ft	1 roll

#### ORTHO-GLASS® Pre-Cut Splints

Order Number	Size	UOM
OG-2PC	2 in x 12 in	10 splints / box
OG-3PC	3 in x 12 in	10 splints / box
OG-335PC	3 in x 35 in	5 splints / box
OG-4PC	4 in x 15 in	10 splints / box
OG-430PC	4 in x 30 in	5 splints / box
OG-5PC	5 in x 30 in	5 splints / box
OG-545PC	5 in x 45 in	5 splints / box
OG-610PC	6 in x 10 in	10 splints / box
OG-6PC	6 in x 30 in	5 splints / box
OG-645PC	6 in x 45 in	5 splints / box





# Portfolio overview

**Actimove® UPPER EXTREMITY SUPPORTS**

Actimove® Shoulder Support	Actimove® Mitella Comfort Arm Sling	Actimove® Elbow Strap	Actimove® Elbow Support	Actimove® Episport® Epicondylitis Clasp
Actimove® Wrist Stabilizer	Actimove® Thumb Stabilizer	Actimove® Alufoam	Actimove® Rhizo Forte Thumb	

**Actimove® LOWER EXTREMITY SUPPORTS**

Actimove® Knee Brace	Actimove® Knee Brace	Actimove® Knee Brace	Actimove® Stabilizer	Actimove® Support
Actimove® GenuMotion Knee Support	Actimove® Post Op Rom	Actimove® Thigh Support	Actimove® PowerMotion Thigh Support	Actimove® PowerMotion Calf Support
Actimove® Standard Walker Low/High	Actimove® Standard Air Walker Low/High	Actimove® Closed Shell Air Walker Low/High	Actimove® Talocast Ankle Brace	Actimove® Ankle Stabilizer
				Actimove® Patella Strap

## Orthopedic Soft Goods

Actimove® is proud to provide medical devices in a broad range of excellent quality products that are reliable, thoroughly tested and tailored to meet the individual patient's needs.

Whether customers doing exercises, recovering from an injury, overcoming daily arthritis or movement-pain related challenges, they can rely on Actimove® to help them with pain relief and provide the support they need.




## Actimove® Shoulder Support

### Extra Pocket for Optional Hot/Cold Pack

- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Extra pain relief when optional hot/cold therapy pack is placed in build-in pocket
- Supports and protects the shoulder while allowing flexibility
- Adjustable chest and arm closure for a custom fit
- Can be worn left or right




	Order Number	Size	Measurements	Color	UOM
	7557310	S	31.5 in. - 41.5 in.	Black	1 ea
	7557311	M	41.5 in. - 51 in.		
	7557312	L	51 in. - 61 in.		

## Actimove® Mitella Comfort Arm Sling

### Professional Line

- Hand positioning through the thumb loop limits the rotation of the forearm while still allowing full finger movement
- Simple to apply and adjust for an individual fit by using the freely rotating O-rings
- Gentle pocket material protects clothes while also protecting casts or bandages from moisture and dirt



	Order Number	Size	Measurements	Color	UOM
	72819-61	S	34 cm - 36 cm	Blue	1 ea
	72819-62	M	37 cm - 43 cm		
	72819-23	L	42 cm - 48 cm		
	72819-64	XL	48+ cm		

## Actimove® Elbow Strap

### Hot/Cold Pack

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Extra pain relief supported by a removable therapeutic hot/cold pack
- Silicone dot lining for secure fit
- Can be worn left or right

	Order Number	Size	Measurements	Color	UOM
	7574230	Universal	11 in.-15 - .75 in.	Black	1 ea



## Actimove® Elbow Support

### Adjustable

- Slip-on design with easy application and adjustment
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Quick and easy to fit with slip-on design and individually adjustable straps
- Can be worn left or right

	Order Number	Size	Measurements	Color	UOM
	7561730	Universal	7.5 in.-14.5 in.	Black	1 ea



## Actimove® EpiSport® Elbow Clasp Adjustable

- Provides support for Epicondylitis, Tennis Elbow, and Muscle Strain
- Dual action compression works simultaneously to treat lateral and medial epicondylitis (forehand and backhand tennis elbow)
- Comfortable shock absorbing liner helps wick moisture away

	Order Number	Size	Measurements	Color	UOM
	7574300	XS	8 in. - 9 in.	Black	1 ea
	7574301	S	9 in. - 10 in.		
	7574302	M	10 in. - 11 in.		
	7574303	L	11 in. - 12 in.		
	7574304	XL	12 in. - 13.5 in.		




## Actimove® Wrist Stabilizer

### Removable Metal Stay, Left or Right

- Anatomically shaped stay for secure stabilization, managing sprains and strains
- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Stabilizes wrist through an anatomically shaped removable stay for left and right hand application
- Quick and easy to fit with slip-on design and individually adjustable straps
- Special comfort in a premium soft material and seamless design
- Can be easily worn under clothing due to low profile
- Can be worn left or right



	Order Number	Size	Measurements	Color	UOM
	7572931	S	5 in. - 6.25 in.	Black	1 ea
	7572932	M	6.25 in. - 7.5 in.	Black	
	7572933	L	7.5 in. - 9 in.	Black	

## Actimove® Alufoam

### Maleable Aluminum Cushioned Finger Strips

- Excellent finger protection
- Lightweight adjustable, trimmable finger splint




Order Number	Size	Measurements	UOM
7295600	S	.5 in. - 18 in.	10/box
7295601	M	.75 in. - 7.5 in.	
7295602	L	1 in. - 9 in.	

## Actimove® Thumb Stabilizer Extra Stays

### Everyday Supports

- Anatomically shaped stays, managing thumb injuries, keeping finger dexterity
- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Stabilization of the thumb through two integrated anatomically shaped stays
- Comfortable to wear because of premium quality soft material
- Firm and adjustable wrist support from powerful elastic strap
- Full finger movement and hand dexterity from special open design



	Order Number	Size	Measurements	Color	UOM
	7563630	S/M	5.5 in. - 7.5 in.	Black	1 ea
	7563631	L/XL	7.5 in. - 9 in.		
	7563632	S/M	5.5 in. - 7.5 in.	Beige	
	7563633	L/XL	7.5 in. - 9 in.		

## Actimove® Rhizo Forte Thumb

### Customizable Thumb Brace

- Low profile, maleable, dirt and water repellant stabilizer for thumb
- Excellent for basal joint arthritis, game keepers and skiers thumb



	Order Number	Size	Hand / Thumb	Color	UOM
	7623800	S			1 ea
	7623801	M	Right	Gray	
	7623802	L			
	7623803	S			1 ea
	7623804	M	Left	Gray	
	7623805	L			

## Actimove® Knee Brace

### Adjustable Horseshoe, Simple Hinges, Condyle Pads

- Hinged knee sleeve for extra stability, adjustable horseshoe for improved patellar tracking
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue and ligament injuries through extra stability with medial and lateral hinges
- Optional patellar support with an adjustable horseshoe
- Stays in place with improved fit through elastic hook and loop support strap
- Easy application with convenient finger-loops
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally
- Can be worn left or right



	Order Number	Size	Measurements	Color	UOM
	7568750	XS	12 in. - 14 in.	Black	1 ea
	7568751	S	14 in. - 16 in.		
	7568752	M	16 in. - 18 in.		
	7568753	L	18 in. - 20 in.		
	7568754	XL	20 in. - 22 in.		
	7568755	XXL	22 in. - 24 in.		
	7568756	XXXL	24 in. - 26 in.		

## Actimove® Knee Brace

### Wrap Around, Polycentric Hinges, Condyle Pads

- Supports easy application for sore knees through front wrap closure, with polycentric hinges designed to match natural movement
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges
- Stays in place with improved fit through elastic hook and loop support straps
- Can be worn left or right



	Order Number	Size	Measurements	Color	UOM
	7311800	XS	12 in. - 14 in.	Black	1 ea
	7311801	S	14 in. - 16 in.		
	7311802	M	16 in. - 18 in.		
	7311803	L	18 in. - 20 in.		
	7311804	XL	20 in. - 22 in.		
	7311805	XXL	22 in. - 24 in.		
	7311806	XXXL	24 in. - 26 in.		

## Actimove® Knee Brace

### Wrap Around, Simple Hinges, Condyle Pads

- Front wrap closure system supports application, simple hinges for extra stability
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges
- Stays in place with improved fit through elastic hook and loop support strap
- Can be worn left or right



	Order Number	Size	Measurements	Color	UOM
	7550010	XS	12 in. - 14 in.	Black	1 ea
	7550011	S	14 in. - 16 in.		
	7550012	M	16 in. - 18 in.		
	7550013	L	18 in. - 20 in.		
	7550014	XL	20 in. - 22 in.		
	7550015	XXL	22 in. - 24 in.		
	7550016	XXXL	24 in. - 26 in.		

## Actimove® Knee Stabilizer

### Adjustable Horseshoe and Stay

- Sleeve design with adjustable horseshoe for improved patellar tracking
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage patellar pain syndrome with an adjustable horseshoe
- Added knee support from medial and lateral spiral stays
- Stays in place with improved fit through elastic hook and loop support strap
- Easy application with convenient finger-loops
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally
- Can be worn left or right



	Order Number	Size	Measurements	Color	UOM
	7245300	XS	12 in. - 14 in.	Black	1 ea
	7245301	S	14 in. - 16 in.		
	7245302	M	16 in. - 18 in.		
	7245303	L	18 in. - 20 in.		
	7245304	XL	20 in. - 22 in.		
	7245305	XXL	22 in. - 24 in.		
	7245306	XXXL	24 in. - 26 in.		



# Actimove® Knee Support

## Closed Patella

- Fast, easy and secure protection for minor knee injuries
- Enjoy sports, as pain relief and healing from contusions, minor Meniscal damage and injuries, swelling from sprains and strains are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Feel secure as knee is protected by a snug fit that stays in place
- Fast & easy slip-on style



Order Number	Size	Measurements	Color	UOM
7558620	XS	12 in. - 14 in.	Navy	1 ea
7558621	S	14 in. - 16 in.		
7558622	M	16 in. - 18 in.		
7558623	L	18 in. - 20 in.		
7558624	XL	20 in. - 22 in.		
7558625	XXL	22 in. - 24 in.	Black	1 ea
7558626	XXXL	24 in. - 26 in.		
7558627	XS	12 in. - 14 in.		
7558628	S	14 in. - 16 in.		
7558629	M	16 in. - 18 in.		
7558630	L	18 in. - 20 in.		
7558631	XL	20 in. - 22 in.		
7558632	XXL	22 in. - 24 in.		
7558633	XXXL	24 in. - 26 in.		



# Actimove® Post Op Rom

## Post Operative Knee Brace

- Telescopic hinge to fit most patients
- Range of motion brace. Set from -10 degrees to 120 degrees

Order Number	Size	Color	UOM
7627855	Universal	Black	1 ea




# Actimove® GenuMotion Knee Support

## Knee Support

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression and a secure fit.
- Contributes to pain relief and healing by reducing edema, through targeted pressure and an intermittent massage effect thanks to the anatomically shaped circular pressure pad.
- For more wearing comfort especially designed soft, smooth Motion-Comfort-Zone that helps to prevent wrinkling in the popliteal area.
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction.
- Added stabilization with great mobility thanks to the lateral spiral stays.
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials.



Order Number	Size	Measurements		Color	UOM
		Above Knee	Below Knee		
7346810	XXS	13 <sup>3</sup> / <sub>4</sub> " - 15" 35 - 38 cm	9 <sup>7</sup> / <sub>8</sub> " - 11" 25 - 28 cm	White	1 ea
7346811	XS	15" - 16 <sup>1</sup> / <sub>8</sub> " 38 - 41 cm	11" - 12 <sup>1</sup> / <sub>4</sub> " 28 - 31 cm		
7346812	S	16 <sup>1</sup> / <sub>8</sub> " - 17 <sup>3</sup> / <sub>8</sub> " 41 - 44 cm	12 <sup>1</sup> / <sub>4</sub> " - 13 <sup>3</sup> / <sub>8</sub> " 31 - 34 cm		
7346813	M	17 <sup>3</sup> / <sub>8</sub> " - 18 <sup>1</sup> / <sub>2</sub> " 44 - 47 cm	13 <sup>3</sup> / <sub>8</sub> " - 14 <sup>5</sup> / <sub>8</sub> " 34 - 37 cm		
7346814	L	18 <sup>1</sup> / <sub>2</sub> " - 19 <sup>3</sup> / <sub>4</sub> " 47 - 50 cm	14 <sup>5</sup> / <sub>8</sub> " - 15 <sup>3</sup> / <sub>4</sub> " 37 - 40 cm		
7346815	XL	19 <sup>3</sup> / <sub>4</sub> " - 20 <sup>7</sup> / <sub>8</sub> " 50 - 53 cm	15 <sup>3</sup> / <sub>4</sub> " - 17" 40 - 43 cm		
7346816	XXL	20 <sup>7</sup> / <sub>8</sub> " - 22" 53 - 56 cm	17" - 18 <sup>1</sup> / <sub>8</sub> " 43 - 46 cm		
7346850	3XL	22" - 23 <sup>1</sup> / <sub>4</sub> " 56 - 59 cm	18 <sup>1</sup> / <sub>8</sub> " - 19 <sup>1</sup> / <sub>4</sub> " 46 - 49 cm		
7346851	4XL	23 <sup>1</sup> / <sub>4</sub> " - 24 <sup>3</sup> / <sub>8</sub> " 59 - 62 cm	19 <sup>1</sup> / <sub>4</sub> " - 20 <sup>1</sup> / <sub>2</sub> " 49 - 52 cm		
7346817	M+	18 <sup>1</sup> / <sub>2</sub> " - 19 <sup>3</sup> / <sub>4</sub> " 47 - 50 cm	13 <sup>3</sup> / <sub>8</sub> " - 14 <sup>5</sup> / <sub>8</sub> " 34 - 37 cm		
7346809	L+	19 <sup>3</sup> / <sub>4</sub> " - 20 <sup>7</sup> / <sub>8</sub> " 50 - 53 cm	14 <sup>5</sup> / <sub>8</sub> " - 15 <sup>3</sup> / <sub>4</sub> " 37 - 40 cm		
7346808	XL+	20 <sup>7</sup> / <sub>8</sub> " - 22" 53 - 56 cm	15 <sup>3</sup> / <sub>4</sub> " - 17" 40 - 43 cm		
7346807	XXL+	22" - 23 <sup>1</sup> / <sub>4</sub> " 56 - 59 cm	17" - 18 <sup>1</sup> / <sub>8</sub> " 43 - 46 cm		
7346852	3XL+	23 <sup>1</sup> / <sub>4</sub> " - 24 <sup>3</sup> / <sub>8</sub> " 59 - 62 cm	18 <sup>1</sup> / <sub>8</sub> " - 19 <sup>1</sup> / <sub>4</sub> " 46 - 49 cm		
7346853	4XL+	24 <sup>3</sup> / <sub>8</sub> " - 25 <sup>5</sup> / <sub>8</sub> " 62 - 65 cm	19 <sup>1</sup> / <sub>4</sub> " - 20 <sup>1</sup> / <sub>2</sub> " 49 - 52 cm		
7346843	XXS	13 <sup>3</sup> / <sub>4</sub> " - 15" 35 - 38 cm	9 <sup>7</sup> / <sub>8</sub> " - 11" 25 - 28 cm		
7346844	XS	15" - 16 <sup>1</sup> / <sub>8</sub> " 38 - 41 cm	11" - 12 <sup>1</sup> / <sub>4</sub> " 28 - 31 cm		
7346845	S	16 <sup>1</sup> / <sub>8</sub> " - 17 <sup>3</sup> / <sub>8</sub> " 41 - 44 cm	12 <sup>1</sup> / <sub>4</sub> " - 13 <sup>3</sup> / <sub>8</sub> " 31 - 34 cm		
7346846	M	17 <sup>3</sup> / <sub>8</sub> " - 18 <sup>1</sup> / <sub>2</sub> " 44 - 47 cm	13 <sup>3</sup> / <sub>8</sub> " - 14 <sup>5</sup> / <sub>8</sub> " 34 - 37 cm		
7346847	L	18 <sup>1</sup> / <sub>2</sub> " - 19 <sup>3</sup> / <sub>4</sub> " 47 - 50 cm	14 <sup>5</sup> / <sub>8</sub> " - 15 <sup>3</sup> / <sub>4</sub> " 37 - 40 cm		
7346848	XL	19 <sup>3</sup> / <sub>4</sub> " - 20 <sup>7</sup> / <sub>8</sub> " 50 - 53 cm	15 <sup>3</sup> / <sub>4</sub> " - 17" 40 - 43 cm		
7346849	XXL	20 <sup>7</sup> / <sub>8</sub> " - 22" 53 - 56 cm	17" - 18 <sup>1</sup> / <sub>8</sub> " 43 - 46 cm		



## Actimove® Thigh Support

### Adjustable

- Adjustable, accommodating fit that stays in place for strains and contusions
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Feel secure as thigh is protected by an accommodating fit that stays in place
- Adjust the level of support for comfort, no matter the activity
- Can be worn left or right




	Order Number	Size	Measurements	Color	UOM
	7560202	Universal	17.75 in. - 31.5 in.	Black	1 ea



## Actimove® PowerMotion Calf Support

- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression.
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle.
- Promotes injured calf muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps.
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone.
- Secure fit due its anatomical shape and adjustable straps.
- Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system.




	Order Number	Size	Measurements	UOM
	7286604	S	11 <sup>1</sup> / <sub>2</sub> " - 13 <sup>3</sup> / <sub>8</sub> " 29 - 34 cm	1 ea
	7286605	M	13 <sup>3</sup> / <sub>8</sub> " - 15 <sup>3</sup> / <sub>4</sub> " 34 - 40 cm	
	7286606	L	15 <sup>3</sup> / <sub>4</sub> " - 18 <sup>1</sup> / <sub>2</sub> " 40 - 47 cm	



## Actimove® PowerMotion Thigh Support

- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression.
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle.
- Promotes injured thigh muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps.
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone.
- Secure fit due its anatomical shape and adjustable straps. Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system.



	Order Number	Size	Measurements	UOM
	7286706	S	16 <sup>1</sup> / <sub>8</sub> " - 17 <sup>3</sup> / <sub>4</sub> " 41 - 45 cm	1 ea
	7286707	M	17 <sup>3</sup> / <sub>4</sub> " - 19 <sup>3</sup> / <sub>4</sub> " 45 - 50 cm	
	7286708	L	19 <sup>3</sup> / <sub>4</sub> " - 22 <sup>1</sup> / <sub>8</sub> " 50 - 56 cm	
	7286709	XL	22 <sup>1</sup> / <sub>8</sub> " - 24 <sup>7</sup> / <sub>8</sub> " 56 - 63 cm	
	7286710	XXL	24 <sup>7</sup> / <sub>8</sub> " - 27 <sup>1</sup> / <sub>2</sub> " 63 - 70 cm	



## Actimove® Standard Walker Low/High

- The Actimove® Standard Walker offer impressive features
- The high-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait



Order Number	Size	Mens Measurements	Womens Measurements	Color	UOM
7627220	XS	under 4.5 in.	under 5.5 in.	Black	1 ea
7627221	S	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627222	M	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627223	L	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627224	XL	12 1/2 in.	13.5 in.		
7627225	XS	under 4.5 in.	under 5.5 in.		
7627226	S	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627227	M	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627228	L	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627229	XL	12.5 in.	13.5 in.		

## Actimove® Standard Air Walker Low/High

- The Actimove® Standard Air Walker offer impressive features
- High-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait



Order Number	Size	Low/High	Mens Measurements	Womens Measurements	Color	UOM
7627230	XS	Low	under 4.5 in.	under 5.5 in.	Black	1 ea
7627231	S	Low	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627232	M	Low	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627233	L	Low	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627234	XL	Low	12 1/2 in.	13.5 in.		
7627235	XS	High	under 4.5 in.	under 5.5 in.		
7627236	S	High	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627237	M	High	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627238	L	High	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627239	XL	High	12.5 in.	13.5 in.		

## Actimove® Closed Shell Air Walker Low/High

- The Actimove® Closed Shell Walker offer impressive features
- High-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait




Order Number	Size	Low/High	Mens Measurements	Womens Measurements	Color	UOM
7627610	XS	Low	under 4.5 in.	under 5.5 in.	Black	1 ea
7627611	S	Low	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627612	M	Low	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627613	L	Low	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627614	XL	Low	12 1/2 in.	13.5 in.		
7627615	XS	High	under 4.5 in.	under 5.5 in.		
7627616	S	High	4.5 in. - 7 in.	5.5 in. - 8.5 in.		
7627617	M	High	7.5 in. - 10.5 in.	8.5 in. - 11.5 in.		
7627618	L	High	10.5 in. - 12.5 in.	11.5 in. - 13.5 in.		
7627619	XL	High	12.5 in.	13.5 in.		

## Actimove® Talocast Ankle Brace

- Engineered for functional performance to help manage pain while providing a comfortable fit and ease of use
- Anatomically shaped stabilizers that support the ankle joint
- Comes with a long lasting and easy to adjust heel strap
- Offers optional cold therapy for extra pain relief


### Actimove Professional Line TaloCast AirGel

	Order Number	Size	Measurements	Color	UOM
	7311900	S/M		Blue / Black	1 ea
	7311901	L/XL			



- Engineered for functional performance to help manage pain while providing a comfortable fit and ease of use
- Anatomically shaped stabilizers that support the ankle joint
- Comes with a long lasting and easy to adjust heel strap

### Actimove Professional Line TaloCast Air


	Order Number	Size	Measurements	Color	UOM
	7311904	S/M		Blue / Black	1 ea
	7311905	L/XL			



## Actimove® Ankle Stabilizer

### Criss-Cross Straps

- Simulating sports taping, increased stability for ankle and arch of foot
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Increased stability around the ankle and arch of foot with criss-cross straps simulating sports taping
- Stays in place during athletic activity with extra support from circular strap
- Can be worn left or right

	Order Number	Size	Color	UOM
	7561130	Universal	Black	1 ea









# Portfolio overview

STRETCHABLE LIGHT-WEIGHT FABRIC DRESSING



Leukoplast® elastic

TRANSPARENT DRESSING WITH ABSORBENT PAD



Leukomed® T Plus

WIDE AREA FIXATION



Leukoplast® Hypafix®

ELASTIC GAUZE BANDAGE



Elastomull®

HYDROPHOBIC MICROBE BINDING FILM DRESSING



Leukomed® Sorbact®

CONFORMING, ADHESIVE FIXATION DRESSING



Cover-Roll® Stretch

WOUND & EYE IRRIGATION



SplashCap®



EyeCap®

## Wound Care



Essity has a full line of wound care products to support your athletes on and off the field.

4.1 STRETCHABLE LIGHTWEIGHT FABRIC DRESSING

## Leukoplast® elastic Fabric Adhesive Bandages

- Absorptive, non-adherent pad
- Secure adhesion for wound protection
- Contours to different areas on the body
- Variety of shapes and sizes



Order Number	Size	Description	UOM
7645617	1/4 in. x 3 in.	Strips 	100/box
7645606	1 in. x 3 in.		100/box
7645604	1 in. x 3 in.		300/box
7645616	7/8 in. in round	Spots 	100/box
7645613	1 1/4 in. in oval		100/box
7645612	1 1/2 in. x 2 in.	Patches 	100/box
7645610	2 in. x 3 in.		50/box
7645605	4 in. x 2 3/4 in.		50/box
7645607	1 1/2 in. x 1 3/8 in.	Knuckles 	100/box
7645618	2 1/2 in. x 1 in.	Toe Shield (non-sterile) 	100/box
7645608	Sm. finger tip	Digits 	100/box
7645609	Lg. finger tip		50/box
7645611	3 in. x 3 in.	4-Wing 	50/box
7646600	—	Eye Occluder REG (non-sterile)	30/box
7646700	—	Eye Occluder JR (non-sterile)	30/box
7253702	—	Cabinet (empty)	—
7253703	—	Cabinet (filled)	—

4.2 TRANSPARENT DRESSING WITH ABSORBENT PAD

## Leukomed® T Plus Transparent, Total-Island Dressing with Absorbent Non-Adherent Pad



- Waterproof and barrier to bacteria
- Skin friendly adhesive allows for reliable dressing retention
- Should clinical conditions allow, the dressing may remain in place for up to 7 days

Order Number	Dressing Size	Wound Pad Size	UOM
7238200	2 in. x 3 in.	1 1/8 in. x 1.5 in.	50/box
7238201	3 1/8 in. x 4 in.	1.5 in. x 2 3/8 in.	50/box
7238202	3 1/8 in. x 6 in.	1.5 in. x 4 3/8 in.	50/box
7238214	4 in. x 7 7/8 in.	2 in. x 6 3/8 in.	50/box
7238203	4 in. x 10 in.	2 in. x 8 1/8 in.	50/box
7238204	4 in. x 11 7/8 in.	2 in. x 9 3/4 in.	50/box
7238205	4 in. x 13 3/4 in.	2 in. x 11 3/4 in.	50/box

4.3 WIDE AREA FIXATION

## Leukoplast® Hypafix® Wide-Area Fixation



- Secure, rapid fixation that is excellent on contoured parts of the body
- Entire wound dressing is protected against contamination
- Convenient cut-to-size roll eliminates unnecessary waste

Order Number	Size	UOM
71443-04	5 cm x 2 m (2 in. x 2 yds)	1 per box/36 boxes
71443-05	10 cm x 2 m (4 in. x 2 yds)	1 per box/12 boxes
71443-06	15 cm x 2 m (6 in. x 2 yds)	1 per box/12 boxes
71443-01	5 cm x 10 m (2 in. x 11 yds)	1 per box/24 boxes
71443-02	10 cm x 10 m (4 in. x 11 yds)	1 per box/24 boxes
71443-03	15 cm x 10 m (6 in. x 11 yds)	1 per box/24 boxes



2.4 ELASTIC GAUZE BANDAGE

## Elastomull®

### Elastic Gauze Bandage

- High density fabric for a soft, rich texture
- Performance and stretchability to help control bleeding and edema
- Precision fit in areas with constant movement
- Helps prevent slippage once bandage is in place



Order Number	Size	UOM / Case
<b>Non-sterile</b>		
7780200	1 in. x 4.1 yds.	24 rolls
7780300	2 in. x 4.1 yds.	12 rolls
7780400	3 in. x 4.1 yds.	12 rolls
7780500	4 in. x 4.1 yds.	12 rolls
<b>Sterile</b>		
7779700	1 in. x 4.1 yds.	12 rolls
7779800	2 in. x 4.1 yds.	12 rolls
7779900	3 in. x 4.1 yds.	12 rolls
7780000	4 in. x 4.1 yds.	12 rolls

4.5 HYDROPHOBIC MICROBE BINDING FILM DRESSING

## Leukomed® Sorbact®

### Wound Protection Dressing with Sorbact® Technology

- Hydrophobic microbe binding film dressing that safely inhibits, binds and removes hydrophobic microbes with a physical mode of action
- Wound contact layer coated with Sorbact® Technology
- Reliable adhesion that is breathable and shower-proof
- Absorbent wound pad
- Should the clinical condition allow, the dressing can be left in place for 7 days
- For use in the management of dry to low exudating wounds, including clean, contaminated, infected wounds, post-operative wounds and trauma wounds



Order Number	Size	UOM
7619914	2 in. x 3 in.	20/box
7619915	3 in. x 4 in.	20/box
7619916	3 in. x 6 in.	20/box
7619917	4 in. x 8 in.	20/box
7619918	4 in. x 10 in.	20/box
7619919	4 in. x 12 in.	20/box
7619920	4 in. x 14 in.	20/box

4.6 CONFORMING, ADHESIVE FIXATION DRESSING

## Cover-Roll® Stretch

### Conforming, Adhesive Fixation Dressing

- Stretchable non-woven polyester material adapts well to body contours
- Transversely, conforms well even to highly contoured parts of the body, helps prevent constriction and congestion
- Secure fixation

Order Number	Size	UOM
45547	2 in. x 2 yds	36 roll/case
45548	4 in. x 2 yds	12 rolls/case
45549	6 in. x 2 yds	12 rolls/case
45550	8 in. x 2 yds	12 rolls/case
45551	12 in. x 2 yds	12 rolls/case
45552	2 in. x 10 yds	12 rolls/case
45553	4 in. x 10 yds	12 rolls/case
45554	6 in. x 10 yds	12 rolls/case
45555	8 in. x 10 yds	12 rolls/case
45556	12 in. x 10 yds	12 rolls/case



4.7 WOUND & EYE IRRIGATION

## SplashCap®

### Irrigate and flush wounds rapidly

- Rapidly flushes contaminants from wounds
- Acute wounds
- Use with standard irrigation bottle not included

Order Number	Color	UOM
10015	Clear	25/case



## EyeCap®

### Irrigate and flush eyes quickly

- Rapidly flushes contaminants from eyes
- Eye contamination
- Use with standard irrigation bottle not included

Order Number	Color	UOM
20025	Clear	25/case





BSN Medical Inc., an Essity company  
5825 Carnegie Blvd., Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325  
To order toll-free: Tel. (+1) 800 537 1063