


**NEW**

# Cutimed® Sorbact® Hydroactive B

For Removal of Necrotic Tissue and Infection Management



## Four Advantages In One Dressing




- Binds bacteria and renders inert
- Hydropolymer gel matrix helps maintain a moist wound environment
- Stimulates autolytic debridement of fibrin and necrotic tissue
- Ideal for absorption of thin exudate

# 4-in-1: Absorb. Hydrate. Debride. + DACC Technology

This innovative hydropolymer gel dressing provides the unparalleled performance of absorbency with balanced hydration. DACC Technology binds bacteria and renders it inert, versus killing like other silver and antimicrobials, which allows for natural wound healing. Combined, these features promote autolytic debridement, reduce the bacterial load, and offer fast, effective wound healing.

## Proof of Wound Healing



Data of clinical observations (n = 100) indicated a noticeable change of phases of wound healing. Inflammatory phase is decreased and the proliferative phase is increased (granulation and epithelialization) after 10 days of Cutimed® Sorbact® Hydroactive treatment.<sup>1</sup>

## Hydration of Peri-Wound Skin

Water and glycerol allow the dressing to hydrate the wound and peri-wound area.



Improvement in peri-wound skin within one week.

## It's All About Balance

Think of the dressing as a cycle of absorbing exudate and stimulating perfusion. Through osmosis, fluid from the wound is pulled into the absorbent polymers. Therefore, this stimulates the wound and draws fluid from surrounding tissues aiding the body's natural process to clear debris and necrotic tissue from the wound.



Cutimed® Sorbact® Hydroactive B Dressing*	Size	Wound Pad Size	REF No.	Dressings/Box	Suggested HCPCS**
	2.8 x 3.3 in. (7 x 8.5 cm)	2 x 2.6 in. (5 x 6.5 cm)	7993300	10	A6245
	5.5 x 5.5 in. (14 x 14 cm)	3.9 x 3.9 in. (10 x 10 cm)	7993302	10	A6245
	7.5 x 7.5 in. (19 x 19 cm)	5.9 x 5.9 in. (15 x 15 cm)	7993303	10	A6246
	5.5 x 9.4 in. (14 x 24 cm)	3.9 x 7.9 in. (10 x 20 cm)	7993304	10	A6246

\*NOTE: Please refer to the product label and/or package insert for full instructions on safe use of these products.

\*\*This document includes codes that might be used to bill for BSN medical products and related services. Each provider will have to verify the appropriate codes for each patient. It is the provider's sole responsibility to determine and submit appropriate codes, charges and modifiers for services rendered. Providers should contact insurers to verify correct coding procedures prior to submitting claims related to any product or service. BSN medical Inc. cannot guarantee coverage or reimbursement with the codes listed in this billing guide. In all cases, providers will need to follow local payer policies for billing and reimbursement.

1 Hallern von, B. (2012, May). Absorbent DACC coated dressings, a new approach in treating wound infection: A multi-centre study. Poster presented at European Wound Management Association (EWMA) Conference, Vienna, Austria.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com  
To order toll-free: BSN medical (+1) 800 552 1157

62061 R1 ©2015 BSN medical Inc. REV 05/15

