


LEUKOTAPE® K

**Permanent massage effect
and pain relief technique
with full range of motion**



THERAPIES. HAND IN HAND.

www.bsnmedical.com

BSN medical

LEUKOTAPE® K

Leukotape® K supports sensory functions.
Wave pattern adhesive allows air to circulate and full mobility is preserved.

Extremely thin and skin-friendly

- 1 Wave pattern adhesive can support the therapeutic effect
- 2 Stretchable, works ideally with the skin's natural elasticity
- 3 Protective liner for better handling
- 4 Polyacrylate adhesive for reliable adhesion
- 5 Backing material of cotton, hard wearing and skin friendly

Patient oriented and patient friendly

- high user comfort
- durable for up to one week's use
- sports and showering present no problem
- easy to remove

Main uses


- headache
- back pain
- lymphedema
- muscular pain

3 sizes and 6 colors

- range to suit all needs and indications
- discreet skin tone or colored

ORDERING INFORMATION

Ref. Code	Length	Width	Bandages p. inner / p. outer	Color
7297810	5.5 yds	1 in	5/50	beige
7297811	5.5 yds	2 in	5/50	beige
7297812	5.5 yds	3 in	5/40	beige
7297815	5.5 yds	1 in	5/50	red
7297816	5.5 yds	2 in	5/50	red
7297817	5.5 yds	3 in	5/40	red
7297820	5.5 yds	1 in	5/50	blue
7297821	5.5 yds	2 in	5/50	blue
7297822	5.5 yds	3 in	5/40	blue
7297823	5.5 yds	2 in	5/50	black
7297824	5.5 yds	2 in	5/50	light blue
7297825	5.5 yds	2 in	5/50	pink



LATEX
FREE
FORMULA

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157