



TAPES AND BANDAGING

FOR ADDED SUPPORT,
PROTECTION AND HEALING

THERAPIES. HAND IN HAND.

BSN medical

an Essity company

LIGHTPLAST® PRO

LIGHT SUPPORT ELASTIC ADHESIVE BANDAGE

Versatile elastic adhesive bandage for all-purpose taping and strapping of ankles, knees, wrists and hands. The unique fabric backing of **Lightplast® Pro** is made of a cotton/spandex blend which helps it to conform to body movement.

PRODUCT BENEFITS


- Consistent unwind and hand tear assures smoother, faster wrapping
- Moisture-resistant rubber adhesive securely holds bandage against sweat and rain
- Air permeable to help prevent skin maceration
- No-waste core allows use of entire roll
- Conforming bandage provides comfortable, light and cool protection

INDICATIONS

- Treatment of sprains and strains
- Prophylactic use against sports injuries
- Retention of dressings

Color	REF No. / BNR No.	Unit of Measure	Size
● White	76952/7695200	48 rolls/case	1 in x 5 yd
	76953/7695300	32 rolls/case	1.5 in x 5 yd
	76954/7695400	24 rolls/case	2 in x 5 yd
	76955/7695500	16 rolls/case	3 in x 5 yd
	76956/7695600	24 rolls/case	2 in x 7.5 yd
	76957/7695700	16 rolls/case	3 in x 7.5 yd
● Black	72880/7288300	24 rolls/case	2 in x 5 yd
	72881/7288100	16 rolls/case	3 in x 5 yd
	72883/7288300	16 rolls/case	3 in x 7.5 yd

Caution: This product contains dry natural rubber.



LIGHTPLAST® COHESIVE

COHESIVE LIGHT SUPPORT BANDAGE

Lightplast® Cohesive is constructed from a combination of polypropylene fibres and elastic yarns and coated with a natural rubber compound to enable it to adhere only to itself. **Lightplast® Cohesive** is versatile and suitable for use where light support or compression is needed.

PRODUCT BENEFITS

- Easy to use: No pins or clips are required to hold the bandage in place
- Unique construction, non-slip coating: Maintains consistent support and compression, bandage will remain in place even during periods of increased movement, e.g., sports
- Strong, light and durable: Can easily be removed and reapplied if necessary
- Flexible substrate: 150–200% stretch allows the bandage to conform to awkward body contours, leading to increased patient comfort and compliance
- Air permeable: Allows the skin to breathe making longer periods of wear more comfortable
- Patient comfort: Adheres only to itself and not to the patient's skin

INDICATIONS

- General support for sprains, strains and soft tissue injuries
- Prophylaxis: Prevention of sports injuries by wrapping targeted areas
- Retention of dressings, splints and hot/cold packs
- Skin which can be damaged or irritated by adhesive bandages: Ideal for use on delicate skin

Color	REF No.	Unit of Measure	Size
● White	7348900	36 rolls/ctn	2 in x 6 yd
	7348901	30 rolls/ctn	3 in x 6 yd
● Black	7348905	36 rolls/ctn	2 in x 6 yd
	7348906	30 rolls/ctn	3 in x 6 yd

Caution: This product contains dry natural rubber.



CO-PLUS® LF

LIGHTWEIGHT COHESIVE BANDAGE

Co-Plus® LF is a flexible cohesive bandage, ideal for use by healthcare professionals and patients with a sensitivity to natural rubber latex.

PRODUCT BENEFITS

- A unique synthetic construction, incorporating a combination of nylon and viscose yarns coated with a latex-free component to provide cohesive properties
- Highly versatile and easy to use – Co-Plus® LF is hand-tearable in both directions, eliminating the need for scissors
- Co-Plus® LF adheres only to itself, not to skin, making the bandage easy to apply and allows effortless pain-free removal
- The unique construction provides consistent, non-slip support
- Co-Plus® LF is able to conform to even the most awkward body contours due to its flexible substrate
- Porous structure which aids skin breathability
- Available in a wide variety of colors, making Co-Plus® LF not only functional, but also attractive to wear

INDICATIONS

- General support and suppression of soft tissue injuries
- Prophylaxis against sports injuries
- Application and retention of dressings, splints, and hot/cold packs, etc.

Color	REF No.	Unit of Measure	Size
● Beige	7210014	60 rolls/pack	1 in x 5 yd
	7210016	36 rolls/pack	2 in x 5 yd
	7210017	24 rolls/pack	3 in x 5 yd
	7210018	18 rolls/pack	4 in x 5 yd
	7210019	12 rolls/pack	6 in x 5 yd
Assorted Colors	7210020	36 rolls/pack	2 in x 5 yd
	7210021	24 rolls/pack	3 in x 5 yd
	7210033	18 rolls/pack	4 in x 5 yd



LEUKOTAPE® P

STRONG, ADHESIVE NON-STRETCH TAPE FOR RIGID SUPPORT

Leukotape® P is a high strength rigid tape with a very strong zinc oxide adhesive, making it ideal for patellofemoral taping. The Leukotape® P Combo Pack is comprised of a roll of Leukotape® P, along with a roll of Cover-Roll® Stretch, Eucerin creme, and adhesive removers.

PRODUCT BENEFITS

- Leukotape® P is specifically designed for the patellofemoral (McConnell) taping technique for correcting patellar position
- Excellent for managing chronic shoulder problems by realigning shoulder-related structures, taping and retraining muscles
- High tensile strength, yet hand-tearable for easy application

INDICATIONS

- Patellofemoral taping
- Rigid joint immobilization and limitation of movement

Color	REF No. / BNR No.	Unit of Measure	Size
● Beige	76168/7616800	30 rolls/case	1.5 in x 15 yd

Caution: This product contains dry natural rubber.



COVER-ROLL® STRETCH

CONFORMING, ADHESIVE FIXATION DRESSING

PRODUCT BENEFITS

- Soft, stretchable non-woven polyester material adapts well to body contours
- Provides light compression without causing constriction or congestion
- The polyacrylate adhesive holds dressing safely and reliably in place

Color	REF No.	Unit of Measure	Size
● White	45547	1 roll/box	2 in x 2 yd
	45548	1 roll/box	4 in x 2 yd
	45549	1 roll/box	6 in x 2 yd
	45550	1 roll/box	8 in x 2 yd
	45551	1 roll/box	12 in x 2 yd
	45552	1 roll/box	2 in x 10 yd
	45553	1 roll/box	4 in x 10 yd
	45554	1 roll/box	6 in x 10 yd
	45555	1 roll/box	8 in x 10 yd
	45556	1 roll/box	12 in x 10 yd



TENSOPLAST®

ELASTIC ADHESIVE BANDAGE

Tensoplast®, formerly Elastoplast®, is an elastic adhesive bandage which is suitable for medium support or compression. Superior stretch and regain provides firm, controlled compression and support. Non-fray edges prevent the bandage from rucking or marking.



PRODUCT BENEFITS

- High adhesive properties
- Provides strong compression and long-term support
- Water repellent
- Air permeable

INDICATIONS

- Treatment of all conditions which require limitation of movement, e.g., sprains, sport injuries
- Prophylactic use in active sport
- Retention of dressings and devices, e.g., post-orthopedic and general surgery or over a muscular wound

Caution: This product contains dry natural rubber.

Tensoplast® Elastic Adhesive Bandage—individually boxed

Color	REF No.	Unit of Measure	Size
● Beige	02598002	36 rolls / pack	1 in × 5 yd
	02599002	36 rolls / pack	2 in × 5 yd
	02600002	36 rolls / pack	3 in × 5 yd
	02601002	36 rolls / pack	4 in × 5 yd
	02602002	12 rolls / pack	6 in × 5 yd
● White	02593002	36 rolls / pack	1 in × 5 yd
	02594002	36 rolls / pack	2 in × 5 yd
	02595002	36 rolls / pack	3 in × 5 yd
	02596002	36 rolls / pack	4 in × 5 yd
	02597002	12 rolls / pack	6 in × 5 yd

Tensoplast® Elastic Adhesive Bandage—boxed packs

Color	REF No.	Unit of Measure	Size
● Beige	04415001	72 rolls / case	2 in × 5 yd
	04416001	48 rolls / case	3 in × 5 yd

Tensoplast® Elastic Adhesive Bandage—boxed packs

Color	REF No.	Unit of Measure	Size
● Beige	04410001	32 rolls / case	1 in × 5 yd
	04412001	24 rolls / case	2 in × 5 yd
	04413001	16 rolls / case	3 in × 5 yd
	04414001	16 rolls / case	4 in × 5 yd

TENSOPLAST® LF

ELASTIC ADHESIVE BANDAGE

Now available latex-free in the quality you trust. Trusted patient comfort and compliance for those sensitive to latex make **Tensoplast® LF** the optimal choice for a wide range of indications.



PRODUCT BENEFITS

- Latex-free fabric
- Permeable material
- Yellow center line

Color	REF No.	Unit of Measure	Size
● White	7623200	12 rolls / box	2 in × 4.9 yd
	7623201	12 rolls / box	3 in × 4.9 yd
	7623202	12 rolls / box	4 in × 4.9 yd

LEUKOTAPE® K

KINESIO THERAPEUTIC ADHESIVE TAPE

Leukotape® K supports sensory functions. Wave pattern adhesive allows air to circulate and full mobility is preserved. It is stretchable, and works ideally with the skin's natural elasticity. It has a protective liner for better handling. **Leukotape® K** has a polyacrylate adhesive for reliable adhesion, while the backing material of cotton is hard wearing and skin friendly.



PRODUCT BENEFITS

- Supports sensory functions
- Wave pattern adhesive allows air to circulate
- Extremely thin and skin-friendly
- Patient-oriented and patient-friendly

INDICATIONS

- Muscular Pain
- Back Pain
- Headache
- Lymphedema

LATEX
FREE
FORMULA



Color	REF No.	Unit of Measure	Size
● Beige	7297810	5 rolls / pack	1 in × 5.5 yd
	7297811	5 rolls / pack	2 in × 5.5 yd
	7297812	5 rolls / pack	3 in × 5.5 yd
● Red	7297815	5 rolls / pack	1 in × 5.5 yd
	7297816	5 rolls / pack	2 in × 5.5 yd
	7297817	5 rolls / pack	3 in × 5.5 yd
● Blue	7297820	5 rolls / pack	1 in × 5.5 yd
	7297821	5 rolls / pack	2 in × 5.5 yd
	7297822	5 rolls / pack	3 in × 5.5 yd
● Black	7297823	5 rolls / pack	2 in × 5.5 yd
● Light Blue	7297824	5 rolls / pack	2 in × 5.5 yd
● Pink	7297825	5 rolls / pack	2 in × 5.5 yd

SUPPORT BANDAGING

DRESSING RETENTION, GENERAL SUPPORT, LIGHT COMPRESSION,
FIXATION OF ORTHO-GLASS® SPLINTS

Ortho-Glass® Color Wraps are excellent for mild compression support. Its controlled stretch retains strength and elasticity and is secured with hook and loop closures. These latex-free bandages are ideal for patients who are sensitive to natural rubber latex. Available in assorted colors and sizes.

**LATEX
FREE**
FORMULA



ORTHO-GLASS® COLOR WRAP

Color	REF No.	Unit of Measure	Size
Color Pack*	CEB-2	18 rolls/pack	2 in x 5 yd
	CEB-3	18 rolls/pack	3 in x 5 yd
	CEB-4	18 rolls/pack	4 in x 5 yd

*Color Pack: gold, green, navy blue, pink, red and royal blue

The **BSN Elastic Bandage** is latex-free, adjustable and constructed of breathable elastic. Ideal for minor sprains and strains. Each bandage is individually wrapped and may be secured by hook and loop closure. Available in a convenient dispenser box.

**LATEX
FREE**
FORMULA



BSN ELASTIC BANDAGE

Color	REF No.	Unit of Measure	Size
● Beige	1037023	10 rolls/box	2 in x 5 yd
	1037033	10 rolls/box	3 in x 5 yd
	1037043	10 rolls/box	4 in x 5 yd
	1037053	10 rolls/box	6 in x 5 yd
	EB-4FLDB	5 rolls/box	4 in x 10 yd
	EB-6FLDB	5 rolls/box	6 in x 10 yd

TO ORDER:

Toll-Free: 800 552 1157 ■ Fax: 800 835 4325

Email: BSN.Orders@BSNmedical.com

BSN medical Inc., an Essity company
5825 Carnegie Blvd., Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325
To order toll-free: Tel. (+1) 800 552 1157



BSNmedical
an Essity company