

Cutimed®

Closing wounds. Together

Sorbion® Sorbact®

CUTIMED® SORBION® SORBACT® SUPER-ABSORBENT DRESSING

FOR REMOVING EXUDATE AND BACTERIA

Superabsorbent core of fibers and gel-forming polymers absorb and retain large volumes of exudate



Wound contact layer coated with Sorbact® Technology safely and effectively binds bacteria without a chemically active agent

Anti-strikethrough backing layer prevents leakage and enables easy dressing changes



THERAPIES. HAND IN HAND.

 **essity**

THE FIRST SUPER-ABSORBENT DRESSING THAT ALSO MANAGES BACTERIA

SAFE REMOVAL OF BACTERIA

- Cutimed® Sorbion® Sorbact® featuring Sorbact® Technology safely and irreversibly binds bacteria, renders them inactive preventing toxin release and replication
- Unlike other antimicrobial wound dressings, bacteria binds to Sorbact® Technology, so cell debris may be safely removed from the wound bed
- An excellent replacement for silver as there is no release of chemically active agents, eliminating the risk of bacterial resistance



ABSORBS AND RETAINS EXUDATE

- The superabsorbent core absorbs and retains moderate to high levels of exudate reducing the risk of skin maceration, even under compression

Ordering Information

| | Size | Ref-No. | Dressings/ Box | Suggested HCPCS** |
|-----------------------------------|------------|---------|-------------------|----------------------|
| Cutimed® Sorbion® Sorbact® | | | | |
| | 4 x 4 in. | 7269808 | 10 | A6196 |
| | 4 x 8 in. | 7269809 | 10 | A6197 |
| | 8 x 8 in. | 7269810 | 10 | A6197 |
| | 8 x 12 in. | 7269811 | 10 | A6198 |


**REIMBURSEMENT DISCLAIMER:

The Noridian Medicare Pricing, Data Analysis and Coding (PDAC) website located at https://www.dmepdac.com/review/items_requiring_coding_verification_reviews.html maintains a searchable database for HCPCS codes based on manufacturer and specific product information. BSN Medical Inc., however, makes no representations as to the accuracy of the information contained within the PDAC database, nor any representations as to whether its products are reimbursable under any government sponsored healthcare program and/or private-insurance program. It is the provider's sole responsibility, in consultation with the insurer, if necessary, to determine and submit appropriate codes, charges, and modifiers for services rendered.

Indications

Cutimed® Sorbion® Sorbact® is indicated for use in the prevention of infection and management of clean, colonized, contaminated or infected, shallow wounds with high to very high exudate levels such as venous leg ulcers, pressure ulcers, diabetic foot ulcers, post-operative and trauma wounds. It can also be used as an absorbent secondary dressing on deeper wounds.

Sorbact® is a registered trademark of ABIGO Medical AB.
Please refer to the product label and/or package insert for full instructions on the safe use of these products.



BSN medical Inc., an Essity company
5825 Carnegie Blvd., Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325
To order toll-free: Tel. (+1) 800 537 1063

63199 R1 ©2018 BSN medical Inc. REV 12/18