


CUTIMED® HYDROCONTROL

FOR REMOVAL OF NECROTIC TISSUE
AND MOIST WOUND HEALING

● Supports Autolytic Debridement

● Optimum Moisture Balance



● Vertical Absorption For Superior Fluid Management

THE PERFECT BALANCE FOR IDEAL MOIST WOUND HEALING

Moisture Management That's Superior In Every Way

- Exceptional fluid management
- Moisture donation
- Ideal for light-to-moderately exuding wounds

Comfort and Patient Satisfaction Elevated


The advanced structure of Cutimed® HydroControl contributes to exceptional comfort, ensuring a high level of patient satisfaction.

- Combination of thinness and high conformability ensures comfort and ease of wear under clothes
- Reliable, skin-friendly adhesive keeps dressing securely in place


Multiple Advantages

Cutimed® HydroControl delivers clinical benefits that save time and reduce costs.

- **Absorption of Exudate** – exudate from the wound bed is absorbed into the hydropolymer gel matrix
- **Movement of Fluid Into the Wound** – fluid moves from the interstitial tissue into the wound due to osmotic effect
- **Autolytic Debridement** – the fresh wound fluid also brings enzymes into the wound that helps to remove fibrin



Excellent results*



Cutimed® HydroControl is able to manage higher levels of fluid **

Cutimed® is a registered trademark of BSN medical GmbH, Comfeel® Plus Transparent is the trademark of Coloplast A/S, Mepilex® Lite is the trademark of Mölnlycke Health Care AB, ALLEVYN® Thin is the trademark of T.J. Smith & Nephew Ltd. and Hydrocoll® Border is the trademark of Paul Hartmann AG. BSN medical GmbH is not associated with any of these companies regarding the trademarks and respective products. Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.

Cutimed® HydroControl

	Size	Ref-No.	Dressings/Box	Shipping Unit	HCPCS
	1 3/4 x 1 3/4 in.	73231-00	10	1/30	A6234
	3 x 3 in.	73231-01	10	1/30	A6234
	4 x 4 in.	73231-02	10	1/10	A6234
	6 x 6 in.	73231-03	10	1/10	A6235

Indications

Cutimed® HydroControl is indicated for dry to moderately exuding, shallow, chronic and secondary healing wounds such as:

- Venous or arterial leg ulcers
- Pressure ulcers
- Diabetic foot ulcers

Cutimed® HydroControl can be used under compression.

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

REIMBURSEMENT DISCLAIMER:

This document includes suggested HCPCS codes that might be used to bill for BSN medical products and related services. Each provider will have to verify the appropriate codes for each patient. These codes are determined by Noridian Healthcare Solutions as the Pricing, Data Analysis and Coding (PDAC) Contractor of the Centers for Medicare and Medicaid Services. These codes may change at the discretion of the PDAC at any time and BSN is not responsible for such changes.

It is the supplier or provider's sole responsibility to determine and submit appropriate codes, charges, and modifiers for services rendered. Providers should contact insurers to verify correct coding procedures prior to submitting claims related to any product or service. BSN medical Inc. cannot guarantee coverage or reimbursement with the codes listed in this billing guide. In all cases, providers will need to follow local payer policies for billing and reimbursement.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.bsnmedical.com

62462 RN ©2016 BSN medical Inc. REV 01/16

