

“ My patients don't want their recovery to hold them back from getting on with life. That's why I use products that are flexible and easy to apply, regardless of wound location. This makes my daily routine easier. ”



FITTING YOUR CARE TO YOUR PATIENT

CUTIMED® SILTEC® – THE BODY-FRIENDLY, ADAPTABLE DRESSING:

Wound dressings have to fit the unique, usually irregular shapes of wounds, conform to the body's contours, and accommodate a certain amount of body movement. At the same time, the dressing has to maintain intimate contact with the wound bed to create an optimal environment for the healing process.

- Comfortable and Conformable**
 The soft and pliable conformability ensures that the dressing can adapt to the large variety of body contours. The body-friendly shapes, with rounded corners, enable a perfect fit wherever the dressing is applied.
- Cutable**
 Cutimed® Siltec® comes in a variety of shapes that fit neatly in any desired position. However, in some circumstances, wound dressings have to be customized to meet special requirements. Cutimed® Siltec® can be cut to any shape without losing its function. The integrated stripes of superabsorbers stay in place and maintain their absorption capacity when the product is cut to size.



Flexibility of Cutimed® Siltec® on the elbow



Flexibility of Cutimed® Siltec® on the heel



Flexibility of Cutimed® Siltec® on the hand

CUTIMED® SILTEC® WOUND CARE THERAPY PRODUCTS

Order information					
Cutimed® Siltec®	Pad size	Wound pad	Dressing size	Ref.-no.	Dressings/box
Cutimed® Siltec® L (Light) Thin silicone foam dressing with super absorbers			2 in x 2.5 in	7328300	10
			4 in x 4 in	7328301	10
			6 in x 6 in	7328302	10
Cutimed® Siltec® Light adhesive foam dressing with super absorbers			2 in x 2.5 in	7328500	10
			4 in x 4 in	7328501	10
			4 in x 8 in	7328502	10
			6 in x 6 in	7328503	10
			8 in x 8 in	7328504	5
Cutimed® Siltec® Plus Gentle and reliable silicone foam dressing with super absorbers			2 in x 2.5 in	7328800	10
			4 in x 4 in	7328801	10
			4 in x 8 in	7328802	10
			6 in x 6 in	7328803	10
			8 in x 8 in*	7328804	5
					*Coming soon
Cutimed® Siltec® B (Border) Fully adhesive silicone border foam dressing with super absorbers	(1.5 x 1.5 in)		3 in x 3 in	7328400	10
	(3 x 3 in)		5 in x 5 in	7328401	10
	(4 x 4 in)		6 in x 6 in	7328402	10
	(5 x 5 in)		7 in x 7 in	7328403	5
	(6.5 x 6.5 in)		9 in x 9 in	7328404	5
	(2.5 x 2.5 in)		4 in x 4 in	7328405	10
	(2.5 x 6.5 in)		4 in x 9 in	7328406	10
The anatomical specialists		Wound pad	Dressing size	Ref.-no.	Dressings/box
Cutimed® Siltec® B Oval Fully adhesive silicone border foam dressing with super absorbers			3 in x 4 in oval	7328414	10
Cutimed® Siltec® Sacrum Fully adhesive silicone border foam dressing for the sacral area with super absorbers			7 in x 7 in	7328700	5
				9 in x 9 in	7328701
Cutimed® Siltec® Plus Heel Fully adhesive silicone foam dressing for the heel with super absorbers			6 in x 9.5 in / Heel 2D	7328805	10
				6 in x 9.5 in / Heel 3D	7328601

Wound Depth	Superficial	Wound Phase	Sloughy	Exudate	Low
	Deep		Granulating		Low to moderate
	Superficial + deep		Epithelialising		Moderate to high

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

BSN medical Inc. 5825 Carnegie Blvd., Charlotte, NC 28209-4633
 Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.bsnmedical.com
 To order toll-free: BSN medical (+1) 800 552 1157

62890 RN ©2017 BSN medical Inc. REV 03/17



Cutimed® Siltec®

Closing wounds. Together.

Super Absorbent Silicone Foam Dressing



What foam dressing excels at helping prevent and manage exudate related challenges?

Cutimed® Siltec® The foam dressing that meets the moisture balance challenge.



THERAPIES. HAND IN HAND.
 www.bsnmedical.com



MANAGING EXUDATE


CUTIMED® SILTEC® TAKES EXUDATE MANAGEMENT TO ANOTHER LEVEL

In order for a wound to heal properly, effective exudate and fluid management is crucial. During the healing process, exudate levels can fluctuate, making it essential that wound dressings provide consistent, reliable results.

DESIGNED TO WORK FOR YOU

THE NEW CUTIMED® SILTEC® COMPONENTS

The right combination of specially designed materials allows Cutimed® Siltec® to deliver the performance you expect from a foam dressing.



CUTIMED® SILTEC® PRODUCTS ENSURE:

- Safe fluid handling**
 Wound exudate has to be absorbed reliably and retained in the foam safely. Cutimed® Siltec® absorbs exudate vertically and transports it to the "hot melt superabsorbent stripes" on top of the foam layer. They expand and form a gel that safely retains a large volume of fluid, even under compression.
- Absorption of viscous exudate**
 Cutimed® Siltec® performs well in handling viscous exudate. This is the result of the perforation of the silicone adhesive layer and the large pores of the foam.
- Protection of the wound edge with minimized risk of maceration**
 Cutimed® Siltec® absorbs fluid vertically, drawing it away from the wound edges, spreading it evenly inside the foam. Additionally, the adhesive silicone layer covers the wound edge and protects it from contact with exudate.




“ When it comes to my patient's care, I choose to use products that I trust and would use on my own loved ones. Their healing is my top priority. ”

TAKING THE PRESSURE OFF

CUTIMED® SILTEC® – A CUSHION WHERE IT IS NEEDED

Immobile patients are often at risk of developing additional skin complications requiring extra attention and care. Cutimed® Siltec® may assist in the prevention of pressure injury as part of a comprehensive plan of continuous care, risk assessment and preventive care by a healthcare facility and healthcare professional.



REDUCING PRESSURE RISK:

- Cushion effect of the foam pad**
 Cutimed® Siltec® helps to distribute pressure and prevents pressure peaks, as shown in the graphic above.
- Low surface friction**
 The smooth surface and rounded edges of the PU film reduce friction between Cutimed® Siltec® and the patient's clothing or bedding.

MANAGING YOUR PATIENT'S RECOVERY

WITH CUTIMED® SILTEC®

The adherent layer of a wound dressing has to hold the dressing securely in place and support the wound healing process. At the same time the dressing has to allow quick and easy handling, so it helps to know that accidental inaccurate application is easy to correct.

CUTIMED® SILTEC® PRODUCTS PROVIDE:



Cutimed® Siltec® silicone is gentle yet provides secure adherence

- Secure fixation**
 Cutimed® Siltec® B (with an adhesive border) has a state-of-the-art silicone layer with an extra secure adhesive formula. It prevents the dressing from shifting and simultaneously helps create an optimal environment for the wound to heal.
- Readjustability**
 While strong enough for secure fixation, the adherent layer in Cutimed® Siltec® can also be removed and reapplied when necessary. This readjustability makes application simple as the dressing can easily and comfortably be repositioned, without wasting dressings.
- Atraumatic removal**
 To ensure that pain at dressing change is minimized, the wound contact layer of Cutimed® Siltec® does not form a physical or chemical bond with the wound bed. When the dressing is removed the wound isn't disturbed and the patient's pain is reduced.

“ It's important I use top-of-the-line products that stay reliably in place but are also easy to adjust and comfortable for my patients. ”