


## CUTIMED® SORBION® SACHET SUPERABSORBENT DRESSINGS FOR WOUND BED PREPARATION

Polypropylene, hypoallergenic outer layer is safe for sensitive skin and maintains structural integrity preventing leakage

Patented design allows space for swelling in absorbent layer that enables optimal dressing contact with the wound bed



Fibers and gelling agents absorb and retain exudate, cell debris and toxins

hydration  
**RESPONSE**  
technology



### Cutimed® Sorbion® Sachet S

- Easy application – either side of dressing may be placed on wound



### Cutimed® Sorbion® Sachet Border

- Self-adhesive border that is flexible, breathable, and waterproof



### Cutimed® Sorbion® Sachet XL

- A shape and size for larger wounds such as lower legs, sacral region, thorax, thigh, and abdomen



### Cutimed® Sorbion® Sachet Multi-Star

- Flexible shape for difficult-to-dress wounds such as toes, heels, sacrum, and elbows



### Cutimed® Sorbion® Sachet S Drainage


- Articulated shape for use with catheters and drainage tubes requiring a slit to ensure a good fit

A variety of shapes and sizes for moderate to very high exuding wounds.

# A SUPERABSORBENT DRESSING WITH MULTIPLE BENEFITS

## OUTSTANDING ABSORPTION CAPACITY

- Securely manages moderate - very high levels of exudate – absorbs 60% more than competitive dressings



**60%**  
more than  
competitive  
dressings






## RETENTION OF EXUDATE

- Retains exudate to reduce maceration, even under compression
- Absorption and retention capacity helps reduce bioburden (trapped within the absorbed exudate) from migrating back into the wound, even under compression

## SOFT DEBRIDMENT

- Dressing creates a moist environment that loosens excess exudate to ensure slough and toxins are gently removed during dressings changes
- Removes the exudate from the wound, which reduces odor

## Ordering Information

	Size	Ref-No.	Dressings/ Box	Suggested HCPCS**
<b>Cutimed® Sorbion® Sachet S</b>				
	3 x 3 in.	7323200	10	A6196
	4 x 4 in.	7323206	10	A6196
	8 x 4 in.	7323209	10	A6197
	8 x 8 in.	7323215	10	A6198
	12 x 8 in.	7323218	10	A6198
<b>Cutimed® Sorbion® Sachet Border</b>				
	4 x 4 in.	7323600	10	A6196
	10 x 6 in.	7323604	10	A6197
<b>Cutimed® Sorbion® Sachet XL</b>				
	18 x 10 in.	7324000	10	A6198
<b>Cutimed® Sorbion® Sachet Multi Star</b>				
	3 in.	7323701	10	A6196
	5.5 in.	7323703	10	A6197
<b>Cutimed® Sorbion® Sachet S Drainage</b>				
	4 x 4 in.	7324100	10	A6196

## Indications

**Cutimed® Sorbion® Sachet** is indicated for moderately to highly exuding wounds such as venous leg ulcers, pressure ulcers, diabetic foot ulcers, dehisced wounds, second-degree burns, skin fistulae, fungating carcinomas, and other types of exuding wounds.

### \*\*REIMBURSEMENT DISCLAIMER:

This document includes suggested HCPCS codes that might be used to bill for BSN medical products and related services. Each provider will have to verify the appropriate codes for each patient. These codes are determined by Noridian Healthcare Solutions as the Pricing, Data Analysis and Coding (PDAC) Contractor of the Centers for Medicare and Medicaid Services. These codes may change at the discretion of the PDAC at any time and BSN is not responsible for such changes.

It is the supplier or provider's sole responsibility to determine and submit appropriate codes, charges, and modifiers for services rendered. Providers should contact insurers to verify correct coding procedures prior to submitting claims related to any product or service. BSN medical Inc. cannot guarantee coverage or reimbursement with the codes listed in this billing guide. In all cases, providers will need to follow local payer policies for billing and reimbursement.

Please refer to the product label and/or package insert for full instructions on the safe use of these products.

All trademarks listed are the property of their respective owners.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325 www.bsnmedical.com