

Delta-Cast[®] Saw & Vacuum

All you ever desired in one compact system



- Convenient and healthy solution
- Superior performance
- Compact all-in-one unit
- Trusted BSN medical service

Based on experience, de



The new **Delta-Cast® Saw** and **Delta-Cast® Vacuum** have been developed with a close insight into the procedures in the cast room. The result of studies and enquiries among cast technicians and practitioners have greatly influenced and inspired the design of this new cast saw and vacuum unit. **Wishes and demands** voiced by key opinion leaders, doctors and cast technicians have been directly incorporated into the design in order to make daily work in the cast room easier, **more convenient and healthy**.

Evaluations after the introduction of **Delta-Cast® Saw** in 2011 show how well the **new design** has been received. Clinicians rated most features (ease of use, weight, noise, vibration, blade change & access) **significantly better** than the saw model they were currently using.

signed with a vision.



A closer look at the features and benefits demonstrates the superior design of Delta-Cast® Saw:

- ✓ Extremely **durable design** and components. Delta-Cast® Saw provides an outstanding lifetime of use.
- ✓ Comfortable, **lightweight design** (1.5 pounds).
- ✓ **Good balance** allows one-handed operation.
- ✓ Special **grip design** enables safe operation even with wet hands or gloves.
- ✓ Two **speed levels** with high oscillation rates provide high cutting performance.
- ✓ **Low noise and vibration** reduces user fatigue and is less disruptive in the cast room.
- ✓ Efficient **dust collection** traps particles that would get airborne.
- ✓ Blade change requires **no removal** of cast saw hood.
- ✓ **2 year warranty** from date of purchase.
- ✓ Vacuum chamber **can be removed** for cleaning as needed.

Delta-Cast® Saw & Vacuum

More than just the sum of its parts

The new Delta-Cast® Saw and Vacuum unit comes with a 2 year warranty. It is a premium choice regarding value for money and longevity. A true and lasting asset in every cast room.

Ordering information

Code	Description	Qty
73481-20	Delta-Cast® Saw (non-vacuum chamber)	1
73481-24	Delta-Cast® Saw (vacuum chamber)	1
73354-20	Delta-Cast® Vacuum	1
73354-00	Delta-Cast® Saw & Vacuum	1
73354-04	Delta-Cast® Vacuum Filter	1
73354-07	Delta-Cast® Vacuum Replacement Fuses	1
73481-10	Vacuum Chamber Replacement Hood	1
73481-13	Replacement Power Supply	1
73481-15	Replacement Intermediate Cord (saw to power supply)	1
73481-01	Saw Replacement Main Cord	1
73481-08	Blade Tool Kit	1
73354-10	Vacuum Main Cord	1

Blades

Code	Description	Size	Qty
73481-30	Power Blade - Polished Steel Plaster, Composites, Plastics	2.5 in.	1
73481-31	Power Blade - PTFE Plaster and Synthetic Cast Materials	2.5 in.	1
73481-32	Power Blade - Titanium Nitride Plaster, Synthetic, Composites or Plastic	2.5 in.	1
73481-33	Power Blade - Dicronite® Plaster, Synthetic, Composites or Plastic	2.5 in.	1

Delta-Cast Vacuum Specifications

Ref. No.	73481-06; 120VAC, 60 Hz 12.3 Amps
Dimensions (Lowered Position)	28.5 in. tall x 21 in. long x 12.5 in. wide
Power Supply Input	110V - 240V, AC Only, 50 - 60 Hz, 120 Watts
Saw Input	9V / 12V DC, 10 Amps
No Load Speed (Low / High)	10,500 / 15,000 OPM*
Duty Cycle	3 / 57 I / O
IP Rating	IPX0

Please refer to the product label and/or package insert for full instructions on the safe use of these products.
www.delta-cast-saw.com • www.delta-cast-vacuum.com

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 www.bsnmedical.com
To order toll-free: BSN medical (+1) 800 552 1157

